
Safe Cut

REF 181130 181140 181145

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'UTILISATION
ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ
ISTRUZIONI PER L'USO
INSTRUÇÕES DE UTILIZAÇÃO
MODO DE EMPLEO
NÁVOD K POUŽITÍ
INŠTRUKCIE NA POUŽÍVANIE


INHALT / CONTENT

Table des matières / Περιεχόμενα / Contenido / Índice/ Contenido / Obsah / Sadržaj:

INHALT / CONTENT 2
DEUTSCH 4
ENGLISH..... 6
FRANÇAIS 8
ΕΛΛΗΝΙΚΑ 10
ITALIANO..... 12
PORTUGUÊS..... 14
ESPAÑOL 16
ČESKY 18
HRVATSKI 20
SLOVENČINA 22
SYMBOLE / SYMBOLS..... 24
INFO 27

ABBILDUNGEN / FIGURES

Figures / Σχήματα / Figure / Imagens / Imágenes / Obrázky / Slike


Abbildung 1 / Figure 1 / Figure 1 / Σχήμα 1 / Figura 1 /
Imagen 1 / Obrázek 1 / Slika 1


Abbildung 2 / Figure 2 / Figure 2 / Σχήμα 2 / Figura 2 /
Imagen 2 / Obrázek 2 / Slika 2


Abbildung 3 / Figure 3 / Figure 3 / Σχήμα 3 / Figura 3 /
Imagen 3 / Obrázek 3 / Slika 3


Abbildung 4 / Figure 4 / Figure 4 / Σχήμα 4 / Figura 4 /
Imagen 4 / Obrázek 4 / Slika 4


Abbildung 5 / Figure 5 / Figure 5 / Σχήμα 5 / Figura 5 /
Imagen 5 / Obrázek 5 / Slika 5


Abbildung 6 / Figure 6 / Figure 6 / Σχήμα 6 / Figura 6 /
Imagen 6 / Obrázek 6 / Slika 6


Abbildung 7 / Figure 7 / Figure 7 / Σχήμα 7 / Figura 7 /
Imagen 7 / Obrázek 7 / Slika 7

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Keep until package is empty

Intended use:

The **Safe Cut Biopsy needle** is used for the puncture and histological extraction of tissue from the bone marrow in the iliac crest region.

Contraindications:

- Haemophilia or similar blood coagulation disorders
- Infection at the site of puncture

Possible known complications:

- Secondary haemorrhage and haematoma
- Infection
- Injury of the skin, soft tissues, nerves or adjacent organs
- Impairment of respiratory functions by the use of analgesics or sedatives

Product description:

The biopsy needle consists of a stylet with a fixed cap (1), a cannula with a grip (2) and a female Luer Lock connector for syringes, a depth stopper (3) to mark the depth of penetration, an insertion aid (4) and an ejector (5) to push out the tissue sample. The ejector has ring markings (6). The additional cap (7) helps to prevent further blood leakage. The **Safe Cut** is equipped with a Lock-Wire (8) for the extraction and securing of the biopsy material.

Warnings:

- The products are **not** intended for **sternal** bone marrow aspiration.
- The products are not suitable for the injection of bone cement.
- For the safe and effective application of bone marrow puncture systems, the physician must have the relevant knowledge, experience and training in the use of this technique on patients.
- In contrast to the standard Jamshidi method, no back-and-forth movements are to be made for the separation of the biopsy material.
- Follow aseptic procedures! Deep wound infections may occur as severe subsequent complications following the application; the removal of which may require a surgical intervention.
- The products should be used in a clinical environment only.
- Sterile disposable gloves are required.
- Caution, risk of injury! Individual components of these products are very sharp and can easily cause injuries.
- These instructions are not intended as a description of the used interventional or surgical techniques. It is the responsibility of the attending physician to ensure the safe use of these products and to consider or evaluate the technique for each individual patient. Do not use if the packaging is damaged or open.
- Never use the product after the expiry date stated on the packaging.
- The needle is individually packaged in transparent rigid blister and medical paper and is sterilized with Ethylene Oxide (EO) in accordance with the requirements of the ISO 11135 standard. It is marketed in packs of 10 pieces. The products are intended explicitly for single use only. The product may not be reused after a single application. The quality of the materials, coatings and adhesive bonds may deteriorate. A safe application is therefore no longer guaranteed. The product is not designed for cleaning and sterilization processes after a single use. The sterility of reprocessed single use products is therefore not guaranteed. The risk of unintentional injuries and infections, especially cross-infection in patients and medical staff considerably increases. The company SOMATEX® Medical Technologies holds no liability in the case of the reesterilisation and/or reuse of the device.
- The products are not designed for MRI safety areas. There is a risk of injury!
- It is possible, that no or too little of the biopsy material is extracted, so that the application has to be repeated.

Preparation:

- The area surrounding the puncture site should be covered to create aseptic conditions.
- Appropriate measures to secure the biopsy sample (biopsy material slides, sample container) for haemopathological evaluation should be taken prior to the intervention.

Biopsy application:

(Follow the illustrations!)

1. Perform subcutaneous and periosteal local anaesthesia at the puncture site.
2. If necessary, use a sharp scalpel to make an incision at the puncture site to facilitate the penetration of the skin.
3. Remove the protective tube from the cannula. The depth stopper (3) should remain on the cannula.
4. **Note:** Prior to the puncture, the depth stopper can be adjusted (3) on the cannula shaft for visual control of the penetration depth (see figure 1).
5. Position the puncture needle (2) and push it forwards until you feel a noticeable resistance indicating that you have reached the bone tissue.
6. Penetrate the needle into the iliac crest in a targeted manner using continuous pressure while rotating (right/left rotation).
7. Upon reaching the medullary canal, the resistance will decrease as it penetrates. Then pull the trocar out of the puncture needle (see figure 2) by lifting off the cover (1).
8. Optionally, blood can now be aspirated from the medullary canal using a syringe.
9. To achieve better handling and to avoid the further leakage of blood, the additional cap (7) should be placed on the needle handle (see figure 2).
10. Now push the needle further into the medullary canal using rotational movements until the required biopsy amount is obtained (see figure 3).
Note: The length of the biopsy sample can be checked by inserting the ejector (5). To do so, remove the additional cap (7) from the handle and insert the ejector (5) carefully into the cannula until you feel noticeable resistance. The cm markings (6) indicate the length of the biopsy sample within the puncture needle (see figure 4).
11. Remove the closing cap (7). Insert the Lock-Wire (8) into the puncture cannula (2) until the handle of the Lock-Wire (8) clicks in (see figure 5).
12. To separate the biopsy from the medullary canal, rotate the complete puncture needle 360° around the cannula shaft at least twice (see Figure 6). The cannula position should not be changed. **Caution: In contrast to the standard Jamshidi method, no back and forth movements are to be made for releasing of the biopsy sample.**
13. Then remove the needle from the bone.
14. Treat the wound.

Collection of the biopsy sample from the needle after the puncture:

1. Remove the Lock-Wire (8) from the puncture needle (1).
2. Place the insertion aid (4) on the tip of the puncture needle in the direction of the arrow.
3. Insert the ejector (5) into the conical opening of the injection aid and carefully remove the biopsy sample from the cannula from the distal to the proximal end (see figure 7).
4. Dispose of the used puncture needle in the appropriate manner.

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SOMATEX Medical Technologies GmbH